

Identification and characterization of small molecules inhibiting RNA targets: miRNA

Julia Frappier¹, Eugenie Laigre², Nathalie Huss-Braun¹, Sophie Darquenne³, Béatrice Vayssière³, Louis Kretz², Clément Soleilhavoup¹, Nathalie Membré¹, Henry Kim², Viviane Brel¹, Renaud Morales²

¹ Department of *in vitro* pharmacology, Novalix, 67000 Strasbourg, France
² Department of biophysics, Novalix, 67000 Strasbourg, France
³ Department of *in vitro* pharmacology, Novalix, 93230 Romainville, France

Abstract

RNA plays a key role in transcription¹ and translation:² gene regulation, proliferation, apoptosis, and development. The biological function is linked to the unique secondary and tertiary structures.³ Many diseases are associated with dysregulation of RNA functions,^{4,5} so there is a growing interest in RNAs as drug targets for novel therapies or targeting undruggable proteins. The most used modalities are antisense oligonucleotides (ASOs), CRISPR gene editing, and more recently small molecules.⁶


Although ASOs and CRISPR are highly specific, they bring some significant challenges such as delivery and adverse reactions. Small molecules are attractive because drug-like properties can be optimized by classical medicinal chemistry.⁷ Targeting RNAs with small molecules is a promising modality for oncology, cardiovascular diseases and CNS disorders. In particular, dys-regulation of microRNAs play key roles in the genesis of the different diseases and thus have been targeted by small molecule inhibitors recently. Here, we describe the setup of

a flowchart to identify and optimize new efficient inhibitors of microRNA. Starting with SPR to select and characterize a group of small molecules that bind specifically to pre-miR with high affinities, we validated their inhibitory activity on miR maturation and function with biochemical and cellular assays. The approach demonstrates a miRNA targeting platform that could be applied to other oncogenic miRNA targets.

Results

HTS screening with SPR: Binding Level Screen

- Annealing
- 3 pre-miRNAs tested in parallel
- Biot-pre-miRNAs immobilization
- Batches of 800-1000 cmpds @ 1 or 2 conc
- Ref compound for set up




Category	Response level (R)
High2, Multibinder	R > 100% Rmax
High	50% Rmax < R < 100% Rmax
Medium	25% Rmax < R < 50% Rmax
Low	R < 25% Rmax
No binding	R ≤ 0

✓ Classification of the binders

Hits characterization and profiling: Multi-Cycle Kinetic

Kinetic fit for Compound A		F02*	
Ch1	miRNA A Seq 1	Ch5	miRNA B Seq 1
Ch2	miRNA A Seq 2	Ch6	miRNA B Seq 2
Ch3	miRNA A Seq 3	Ch7	miRNA B Seq 3
Ch4	miRNA A Seq 4	Ch8	miRNA B Seq 4

- Affinity & Kinetic determination
- Hits ranking K_D
- Hits profiling: selectivity vs other miRNAs



High-throughput SPR: simultaneous analysis of multiple targets and compounds

	One channel-based SPR	High-throughput SPR
Time (50 compounds run in dose-response)	70 hrs	15 hrs
Compounds per assay	Max 384/ 1 cmpd per cycle	32Pro = 12 x 384 8K = 4 x 384 8 cmpds per cycle in parallel
Number of targets analysed in parallel	3	32Pro: up to 24 8K: up to 8

Specificity assessment by SPR


- Counter screen performed against Ribonuclease of interest
- Kinetic analysis reveals subtle differences in binding between DICER and miRNA providing further insights into MoA

1 - Fluorescent assay to measure DICER activity → 2 - Validation of the miRNA cleavage by Bioanalyzer → 3 - Cellular activity (gene reporter): target engagement


Fluorescent quenching assay principle



Assay development in an HTS format for 2 pre-miRNA (mutant vs wt)



Compound screening at 2 doses then DRC for IC50 determination




Fluorescence quenching assay


- PROS**
- High to medium throughput
 - IC50 determination
 - Low costs
- CONS**
- Fluorescent cmpds interference
 - No information about cleavage specificity

Analysis of cleaved miRNA

Compound screening LTS: 4 chips (11 channels/day)
 Assay compatible for autofluorescent compounds



Quantification




Pre-miRNA cleavage w/ Bioanalyzer


- PROS**
- Quantification of % of cleavage
 - No interference with fluorescent cmpds
 - Cleavage specificity: RNA pieces size giving insights about MoA
- CONS**
- Low throughput
 - Expensive

psiCHECK2 gene reporter assay development on 5 cell lines for 2 different miRNA

HTS format




Compound screening at 2 doses then DRC for IC50 determination




4 - Antiproliferative activity & toxicity

Cell Titer Glow proliferation assay (panel of cancer cells)




Toxicity is not linked to miRNA inhibition

Toxicity on normal cells



Off-target toxicity elucidation is ongoing

Flowchart



Summary

An integrated workflow for the identification of small molecule inhibitors of RNA targets is demonstrated. High-throughput screening by SPR allowed rapid identification of hit compounds and the characterization of their kinetics and affinity. Biochemical assays showed specific inhibitory activity of pre-miR maturation. Cellular assays confirmed cellular penetration, on-target activity, and anti-proliferation. The workflow was optimized to go rapidly from the hit identification to lead optimization. Similar workflows have been developed for other RNA targets including aptamers, messenger-RNA, ribosomes, RNA-protein complexes, and trinucleotide expansion.

References

- Zamfi B., Bintu L., Ishibashi T. and Bustamante C. (2012) *Proc. Natl. Acad. Sci. U.S.A.* 109, 8948–8953;
- Mustoe A.M., Corley M., Laederach A. and Weeks K.M. (2018) *Biochemistry* 57, 3537–3539;
- Vicens, Q., and Kieft, J.S. (2022). *Proc. Natl. Acad. Sci. USA* 119;
- Bernat V. and Disney M.D. (2015) *Neuron* 87, 28–46;
- Wang F., Zuroski T. and Watts J.K. (2020) *Nat. Rev. Drug Discov.* 19, 441–442;
- Crooke, S. T., Baker, B. F., Crooke, R. M. & Liang, X. H. *Nat. Rev. Drug Discov.* 20, 427–453 (2021);
- Duca, M. et al. *Nat Rev Chem* 8, 120–135 (2024).