

Tailored Fragment Screening Strategies for Diverse Fragments



Albane Marck, Anna Bonhoure, Camille Kostmann, Elise Del Nero, Eugenie Laigre, Johann Stojko, Laura Carneiro, Louis Kretz, Francois Debaene, Henry Kim, Marc Vitorino, Renaud Morales

Department of biophysics, Novalix, 67000 Strasbourg, France

Abstract Fragment-based drug discovery (FBDD) identifies low MW compounds that bind weakly to druggable protein targets. The key hit-finding techniques are surface plasmon resonance (SPR) and nuclear magnetic resonance (NMR). In contrast, covalent binders (with electrophilic warheads) are most found using MS combining intact mass by RPLC-MS (overall distribution of bound ligands), and peptide mapping (site distribution/specificity).

Efficient screening processes rapidly generate hits, that are confirmed by more detailed measurements and orthogonal techniques. Examples of each are presented

Non covalent Fragment screening



Fragment Library

External or internal library

Example: NovaFrag library

- 1200 cpds
- Fully qualified library
- Good coverage of chemical space
- Significant 3D/sp3
- Synthetically tractable
- Rapid follow-up by SAR by catalog
- IP free

Screening Capacity

SPR Cytiva Biacore 8K:

- 1 week screen: singletons at 1 conc
- 1 week screen: hits at 3 concs
- 1 week dose-response: K_D

NMR Bruker Ultra High-field 600 MHz –1 GHz:

- Cocktails of 5 cpds:
- Ligand observed
- NMR experiments: STD, WaterLOGSY, T2/T1rho
- Fully automated: SampleJet, TECAN robot
- 2 weeks for set-up
- Min 4 weeks for screening



Target Carboxylases







- Parallel LO-NMR (cocktails) & SPR screen
- Client library of 2000 fragments
- 700+ NMR detected fragment
- Hit rate ~35 % (unusual, high)
- 20 common HITS between SPR and NMR

Target Acetyltransferase

- Parallel LO-NMR & SPR screen
- 1250 NovaFrag library
- 31 common detected fragment HITS
- Two successive competition assays by NMR
- Kinetic profiling of selected hits by SPR
- 3 categories defined

Target Phosphatase

- 1D ¹H-NMR, Protein Observed
- Mixture of 12 Fragments
- Hits validated with ¹³C-CYS-labeled protein



Covalent Fragment screening

LC-MS: Intact mass




Incubation optimization

- Target residue accessibility
- Warhead reactivity (chloroacetamides, acrylamides etc.)
- Protein stability

PL1 Binding (%)											
ACR	Mean	CIaL	ACR	Mean	CIaL	ACR	Mean	CIaL	ACR	Mean	CIaL
SPR_10 μM	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
SPR_10 μM	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
SPR_100 μM	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%


Throughput ~250 cpds/day

- Client or commercial library
- Automated data processing: UniDec + KNIME
- Ranking of compounds
- Counterscreen on mutated or subdomain of the protein




Binding site identification

- Bottom-up approach
- Quantification
- Determination of site specificity



k_{inact}/K_i

- MedChem follow-up



$$PL1 \text{ binding (\%)} = \frac{I(PL1)}{\sum I(P+PLn)}$$

$$PLn \text{ binding (\%)} = \frac{\sum I(PLn)}{\sum I(P+PLn)}$$

$$DTR = \frac{\sum I(PLn \times n)}{\sum I(PLn)}$$

- Grey: (PL1 = PLn = DTR = 0) and Orange: (PLn ≤ 15%, DTR < 0.3) or (PLn = 100%, DTR > 2) → discarded
- R1 (Rank 1): (PL1 ~ PLn > 50% and 0.7 ≤ DTR ≤ 1.5) → Rank 1 cpds
- R2 (Rank 2): (PL1 ~ PLn > 25% and 0.3 ≤ DTR ≤ 0.6) → Rank 2 cpds

Summary

The examples illustrate efficient fragment screening workflows based on several biophysical techniques. Non-covalent fragment hits were found by screening a chemically diverse library. SPR and NMR each generate some unique hits but may also be used as orthogonal means of hit confirmation prior to SAR enablement by X-ray crystallography and hit-to-lead programs. Covalent fragment screening by MS has been applied to multiple targets with a wide variety of warheads leading to the identification and characterization of specific and potent covalent binders.